

# Laparotomy Sponge (prewashed with loop)



## **Instructions For Use:**

- Offered to O.R team for absorbing blood and exudates or for cleaning wound.
- Prevent the wound from being infected.

## **Features:**

- Material: 100% cotton absorbent gauze.
- 4ply, 6ply, 8ply or 12ply.
- Size: 30x30cm, 45x45cm, 40x45cm, 20x30cm, 20x10cm, 35x35cm, 25x25cm, 30x35cm.
- Prewashed material, with loop.
- With or without x ray.
- Sterile or unsterile.
- Pure white, soft and high absorbency.

Size	Unsterile/Ref	Sterile/Ref	Material
30x30cm, 4ply	190-N1-1-16B	190-S1-1-16B	cotton absorbent gauze
45x45cm, 4ply	190-N1-2-16B	190-S1-2-16B	cotton absorbent gauze
30x30cm, 6ply	190-N2-1-16B	190-S2-1-16B	cotton absorbent gauze
45x45cm, 6ply	190-N2-2-16B	190-S2-2-16B	cotton absorbent gauze
30x30cm, 8ply	190-N3-1-16B	190-S3-1-16B	cotton absorbent gauze
45x45cm, 8ply	190-N3-2-16B	190-S3-2-16B	cotton absorbent gauze
30x30cm, 12ply	190-N4-1-16B	190-S4-1-16B	cotton absorbent gauze
45x45cm, 12ply	190-N4-2-16B	190-S4-2-16B	cotton absorbent gauze
30x30cm, 4ply, with X ray	190-NX1-1-16B	190-SX1-1-16B	cotton absorbent gauze
45x45cm, 4ply, with X ray	190-NX1-2-16B	190-SX1-2-16B	cotton absorbent gauze
30x30cm, 6ply, with X ray	190-NX2-1-16B	190-SX2-1-16B	cotton absorbent gauze
45x45cm, 6ply, with X ray	190-NX2-2-16B	190-SX2-2-16B	cotton absorbent gauze
30x30cm, 8ply, with X ray	190-NX3-1-16B	190-SX3-1-16B	cotton absorbent gauze
45x45cm, 8ply, with X ray	190-NX3-2-16B	190-SX3-2-16B	cotton absorbent gauze
30x30cm, 12ply, with X ray	190-NX4-1-16B	190-SX4-1-16B	cotton absorbent gauze
45x45cm, 12ply, with X ray	190-NX4-2-16B	190-SX4-2-16B	cotton absorbent gauze



**Caution:**

1. This product is sterilized for single use only, discard after single use.

Please do not repeat sterilization, do not use it once again.

2. Should be used the product immediately after opening the individual package, please destroy when it was used.

3. Do not use it if the inner package has been previously opened or damaged, or package be affected with damp, or the product went moldy.

4. Do not use after expiry date. If over the expiry date, please discard it.

5. Sterile by gamma, or non-sterile. Single use only. Latex free.

**Storage:**

Store in a cool, dry and dark place. Avoid direct sunlight.

# Laparotomy Sponge (prewashed with loop)



**Secured Medical Direction UK Co., Ltd.**

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



**EC REP** Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

[www.smd-medical.com](http://www.smd-medical.com)