

Non-woven Swab



Instructions For Use:

- Ideal for wound cleaning and changing for the patient who suffer trauma.
- Absorbing the blood and exudates, haemostasis and wound isolation.
- X-ray tread used for detecting.
- Used in home, hospital, health care, beauty care and so on.

Features:

- Material: non-woven.
- 4ply, 6ply or 8ply.
- Size: 5x5cm, 7.5x7.5cm, 10x10cm, 15x15cm, 10x20cm.
- With or without X ray.
- Sterile or unsterile.

Size	Unsterile/Ref	Sterile/Ref	Material
5.0cm x 5.0cm, 4ply	80-N4-1-16B	80-S4-1-16B	non-woven
7.5cm x 7.5cm, 4ply	80-N4-2-16B	80-S4-2-16B	non-woven
10.0cm x 10.0cm, 4ply	80-N4-3-16B	80-S4-3-16B	non-woven
15.0cm x 15.0cm, 4ply	80-N4-4-16B	80-S4-4-16B	non-woven
10cm x 20cm, 4ply	80-N4-5-16B	80-S4-5-16B	non-woven
5.0cm x 5.0cm, 4ply, with X ray	80-NX4-1-16B	80-SX4-1-16B	non-woven
7.5cm x 7.5cm, 4ply, with X ray	80-NX4-2-16B	80-SX4-2-16B	non-woven
10.0cm x 10.0cm, 4ply, with X ray	80-NX4-3-16B	80-SX4-3-16B	non-woven
15.0cm x 15.0cm, 4ply, with X ray	80-NX4-4-16B	80-SX4-4-16B	non-woven
10cm x 20cm, 4ply, with X ray	80-NX4-5-16B	80-SX4-5-16B	non-woven
5.0cm x 5.0cm, 6ply	80-N6-1-16B	80-S6-1-16B	non-woven
7.5cm x 7.5cm, 6ply	80-N6-2-16B	80-S6-2-16B	non-woven
10.0cm x 10.0cm, 6ply	80-N6-3-16B	80-S6-3-16B	non-woven
15.0cm x 15.0cm, 6ply	80-N6-4-16B	80-S6-4-16B	non-woven
10cm x 20cm, 6ply	80-N6-5-16B	80-S6-5-16B	non-woven
5.0cm x 5.0cm, 6ply, with X ray	80-NX6-1-16B	80-SX6-1-16B	non-woven
7.5cm x 7.5cm, 6ply, with X ray	80-NX6-2-16B	80-SX6-2-16B	non-woven
10.0cm x 10.0cm, 6ply, with X ray	80-NX6-3-16B	80-SX6-3-16B	non-woven
15.0cm x 15.0cm, 6ply, with X ray	80-NX6-4-16B	80-SX6-4-16B	non-woven
10cm x 20cm, 6ply, with X ray	80-NX6-5-16B	80-SX6-5-16B	non-woven
5.0cm x 5.0cm, 8ply	80-N8-1-16B	80-S8-1-16B	non-woven
7.5cm x 7.5cm, 8ply	80-N8-2-16B	80-S8-2-16B	non-woven
10.0cm x 10.0cm, 8ply	80-N8-3-16B	80-S8-3-16B	non-woven
15.0cm x 15.0cm, 8ply	80-N8-4-16B	80-S8-4-16B	non-woven
10cm x 20cm, 8ply	80-N8-5-16B	80-S8-5-16B	non-woven
5.0cm x 5.0cm, 8ply, with X ray	80-NX8-1-16B	80-SX8-1-16B	non-woven
7.5cm x 7.5cm, 8ply, with X ray	80-NX8-2-16B	80-SX8-2-16B	non-woven
10.0cm x 10.0cm, 8ply, with X ray	80-NX8-3-16B	80-SX8-3-16B	non-woven
15.0cm x 15.0cm, 8ply, with X ray	80-NX8-4-16B	80-SX8-4-16B	non-woven
10cm x 20cm, 8ply, with X ray	80-NX8-5-16B	80-SX8-5-16B	non-woven

**Caution:**

1. This product is sterilized for single use only, discard after single use.
Please do not repeat sterilization, do not use it once again.
2. Should be used the product immediately after opening the individual package,
please destroy when it was used.
3. Do not use it if the inner package has been previously opened or damaged,
or package be affected with damp, or the product went moldy.
4. Do not use after expiry date. If over the expiry date, please discard it.
5. Sterile by ETO, or non-sterile. Single use only. Latex free.

Storage:

Store in a cool, dry and dark place. Avoid direct sunlight.

Non-woven Swab



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



EC REP Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com