

# Disposable Dressing Pack

**Disposable Dressing Pack  
(standard pack 1)**



**Disposable Dressing Pack  
(standard pack 2)**



**Disposable Dressing Pack  
(standard pack 3)**



**Disposable Dressing Pack  
(pack L)**



**Disposable Dressing Pack  
(pack M)**



**Disposable Dressing Pack  
(pack M with scissor)**



## Instructions For Use:

- It consists of different accessories.
- All are grouped and placed in a unique pack.
- Idea for wound cleaning and changing for the patients who suffer trauma.
- Used for different department in hospital or home care.

### Name of parts for Disposable Dressing Pack

Disposable Dressing Pack		
Item No.	Name	Material
1	Hard tray	PP
2	Sterile field	PE
3	Patient Drape	Wood pulp, PE
4	Limpet bag	PE
5	Dressing forceps	PS
6	Gauze swab	Cotton
7	Non-woven swab	Non-woven
8	Cotton balls	Cotton
9	Hand towel	Wood pulp
10	Absorbent cotton pad	Non-woven, cotton, wood pulp
11	Plastic handle scissor	Plastic with stainless steel

### Disposable Dressing Pack (standard pack 1)



### REF: 96-1-16A

- 1pc hard tray with 3 compartments.
- 1pc sterile field.
- 1pc patient drape.
- 1pc limpet bag.
- 1pc hand towel.
- 3pcs dressing forceps.
- 4pcs gauze swabs 8ply 7.5x7.5cm.
- 6pcs cotton balls.


**Disposable Dressing Pack  
(standard pack 2)**



**REF: 96-2-16A**

- 1pc hard tray with 3 compartments.
- 1pc sterile field.
- 1pc patient drape.
- 1pc limpet bag.
- 1pc hand towel.
- 1pair latex glove.
- 3pcs dressing forceps.
- 4pcs gauze swabs 8ply 7.5x7.5cm.
- 6pcs cotton balls.

**Disposable Dressing Pack  
(standard pack 3)**



**REF: 96-3-16A**

- 1pc hard tray with 3 compartments.
- 1pc sterile field.
- 1pc patient drape.
- 1pc limpet bag.
- 4pcs dressing forceps.
- 5pcs gauze swabs 8ply 7.5x7.5cm.
- 6pcs cotton ball.
- 2pcs hand towels.

**Disposable Dressing Pack  
(pack L)**



**REF: 96-L-16**

- 1pc hard tray with 3 compartments.
- 1pc sterile field.
- 1pc patient drape.
- 1pc limpet bag.
- 4pcs dressing forceps.
- 5pcs non-woven swabs 4ply 10.5x10cm.
- 10pcs cotton ball.
- 2pcs hand towels.
- 1pc absorbent pad.

**Disposable Dressing Pack  
(pack M)**



**REF: 96-M-16**

- 1pc hard tray with 3 compartments.
- 1pc sterile field.
- 1pc patient drape.
- 1pc limpet bag.
- 2pcs dressing forceps.
- 5pcs non-woven swabs 4ply 10.5x10cm.
- 6pcs cotton balls.
- 2pcs hand towels.
- 1pc absorbent pad.

**Disposable Dressing Pack  
(pack M with scissor)**



**REF: 96-M-S-16**

- 1pc basic dressing tray with 3 compartments.
- 2pcs dressing forceps.
- 1pc limpet bag with tie.
- 6pcs cotton balls.
- 1pc patient drape.
- 5pcs non-woven swabs 4ply 10.5x10cm.
- 2pcs hand towels.
- 1pc absorbent pad.
- 1pc sterile field.
- 1pc plastic handle scissor.

**Disposable Dressing Pack  
(pack S)**



**REF: 96-S-16**

- 1pc basic dressing tray with 3 compartments.
- 1pc dressing forceps.
- 1pc limpet bag with tie.
- 3pcs cotton balls.
- 1pc patient drape.
- 5pcs non-woven swabs 4ply 7x7cm.
- 1pc hand towel .
- 1pc absorbent pad.
- 1pc sterile field.

**Caution:**

1. This product is sterilized for single use only, discard after single use.

Please do not repeat sterilization, do not use it once again.

2. Should be used the product immediately after opening the individual package, please destroy when it was used.

3. Do not use it if the inner package has been previously opened or damaged, or package be affected with damp, or the product went moldy.

4. Do not use after expiry date. If over the expiry date, please discard it.

5. Sterile by ETO. Single use only.

**Storage:**

Store in a cool, dry and dark place. Avoid direct sunlight.

# Disposable Dressing Pack



**Secured Medical Direction UK Co., Ltd.**

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.

**EC REP** Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

[www.smd-medical.com](http://www.smd-medical.com)