

Self-Sealing Sterilization Pouch



Instructions For Use:

- Apply to hospital, clinic and laboratory's sterilization, and also apply to the disinfection of beauty products or family high-temperature.
- Apply to both ETO sterilization and steam sterilization.
- The sharp part of equipment should be put contrary to the peel side to ensure safety use.
- Tear off the release paper aside the tangent line, put the equipment into the pouch, seal the pouch.
- Mark the sterilization time, contents and usage etc.

Features:

- Material: high temperature dialysis paper/France Arjo Paper + medical CPP/PET compound film + medical adhesive tape.
- The medical adhesive tape can make rapid and efficient seal without any professional sealing devices.
- The indicator's color has significant change and no change back. The initial color will change to two different color after steam sterilization and ETO sterilization.
- The three sided explosion-proof can efficiently prevent the pouch from bursting.
- The transparent compound film can show what inside the pouch clearly in order to reduce errors.

Type	Size	Ref	Material
Self-sealing	57x130mm	328-1-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	70x260mm	328-2-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	83x160mm	328-3-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	90x135mm	328-4-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	90x170mm	328-5-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	90x260mm	328-6-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	135x260mm	328-7-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	135x280mm	328-8-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	135x290mm	328-9-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	190x360mm	328-10-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	250x370mm	328-11-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	250x400mm	328-12-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	260x410mm	328-13-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	300x410mm	328-14-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	300x460mm	328-15-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	305x430mm	328-16-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape
Self-sealing	190x330mm	328-17-20B	dialysis paper/France Arjo Paper, CPP/PET compound film, adhesive tape



Caution:

1. Read instruction for use before using the device.
2. This product is for single use only, discard after single use. Please do not use it once again, please destroy when it was used.
3. Do not use it if the package be affected with damp, or the product went moldy.
4. Do not use after expiry date. If over the expiry date, please discard it.
5. Non-sterile. Single use only. Latex free.

Storage:

1. Store in a cool, dry and dark place. Avoid direct sunlight.
2. The clean area with temperature below 25°C and humidity below 60% is recommended.

Self-Sealing Sterilization Pouch



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



EC REP Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com