

Heat Sealing Sterilization Gusseted Pouch



Instructions For Use:

- Apply to hospital, clinic and laboratory's sterilization, and also apply to the disinfection of beauty products or family high-temperature.
- Apply to both ETO sterilization and steam sterilization.
- The sharp part of equipment should be put contrary to the peel side to ensure safety use.
- Using sealing machine to seal and mark the sterilization time, contents and usage etc.

Features:

- Material: imported medical paper coating with high temperature glue.
- The pouch can be sealed by thermal sealing machines.
- Gusseted Pouch, to increase increase the volume of the pouch.
- Printed with sterilization indicators to check the sterilized result.
- Apply to various medical materials and devices.
- Size can be as customer's request.

Size	Ref	Material
230x350x100mm	372-1-20B	medical paper, high temperature glue
250x460x50mm	372-2-20B	medical paper, high temperature glue
100x230x50mm	372-3-20B	medical paper, high temperature glue



Caution:

1. Read instruction for use before using the device.
2. This product is for single use only, discard after single use. Please do not use it once again, please destroy when it was used.
3. Do not use it if the package be affected with damp, or the product went moldy.
4. Do not use after expiry date. If over the expiry date, please discard it.
5. Non-sterile. Single use only. Latex free.

Storage:

1. Store in a cool, dry and dark place. Avoid direct sunlight.
2. The clean area with temperature below 25°C and humidity below 60% is recommended.

Heat Sealing Sterilization Gusseted Pouch



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com