


PE Apron




Instructions For Use:

- Peel off the bag along with dotted line.
- Take out the PE apron.
- Wear on the neck, and tie the waist back.
- For isolation and protective use in different field.

Features:

- Material: polyethylene (LDPE or HDPE).
- Size: 48x28 inch, 53.5x28 inch.
- Style: embossed or smooth.
- Color: white, blue or others.



Size	Ref	Material
L48xW28 inch	137-1-20A	Polyethylene
L48xW28 inch	137-1-20B	Polyethylene
L48xW28 inch	137-1-20C	Polyethylene
L53.5xW28 inch	137-2-20C	Polyethylene

Caution:

1. Read instruction for use before using the device.
2. This product is for single use only, discard after single use. Please do not clean, do not use it once again. Please destroy when it was used.
3. Do not use it if the package be affected with damp, or the product went moldy.
4. Do not use after expiry date. If over the expiry date, please discard it.
5. Non-sterile. Single use only. Latex free.

Storage:

Store in a cool, dry and dark place. Avoid direct sunlight.

PE Apron



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com