

Medical Bedsheet Roll



Instructions For Use:

- Peel off the packing, take out the medical bedsheet roll.
- Fix the bedsheet roll on the bed.
- Tear off piece by piece from the pre-cut place. Cover on the bed.
- Designed for hygiene, cleaning, and protection purpose, avoid cross-infection.

Features:

- Material: wood pulp paper.
- 50cmx99m, 33cm pre-cut.
- 2 layer, embossed.
- Can be 1 layer, 2 layer or 3 layer.
- With PE layer or without.
- Different size can be as required.

Size	Ref	Material
50cmx99M+33cm pre-cut	392-1-20A	Wood pulp paper



Caution:

1. Read instruction for use before using the device.
2. This product is for single use only, discard after single use. Please do not use it once again.
Please destroy when it was used.
3. Do not use it if the package be affected with damp, or the product went moldy.
4. Do not use after expiry date. If over the expiry date, please discard it.
5. Non-sterile. Single use only. Latex free.

Storage:

Store in a cool, dry and dark place. Avoid direct sunlight.

Medical Bedsheet Roll



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com