


Silk Tape



Instructions For Use:

- A strong tape for affixing. Zigzag edges facilitate tearing.
- It is multipurpose use as fixing IV cannula, catheters and other wound dressing.
- Cut the size as per request. Fixed on the place where needed.

Features:

- Material: Rayon silk or acetate silk.
- Zigzag edge.
- Size: 1.25cm, 2.5cm, 5cm, 7.5cm.
- Length: 4.6m, 9.2m.
- Reliable adhesiveness, low sensitization, excellent compliance, no residue glue.

Size	B/Ref	C/Ref	Material
12.5mm(1/2") x10y/9.2m	67-9.2-1-16B	67-9.2-1-16C	Rayon silk or acetate silk
25.0mm(1") x 10y/9.2m	67-9.2-2-16B	67-9.2-2-16C	Rayon silk or acetate silk
50.0mm(2") x 10y/9.2m	67-9.2-3-16B	67-9.2-3-16C	Rayon silk or acetate silk
75.0mm(3") x 10y/9.2m	67-9.2-4-16B	67-9.2-4-16C	Rayon silk or acetate silk
12.5mm(1/2") x5y/4.6m	67-4.6-1-16B	67-4.6-1-16C	Rayon silk or acetate silk
25.0mm(1") x 5y/4.6m	67-4.6-2-16B	67-4.6-2-16C	Rayon silk or acetate silk
50.0mm(2") x 5y/4.6m	67-4.6-3-16B	67-4.6-3-16C	Rayon silk or acetate silk
75.0mm(3") x 5y/4.6m	67-4.6-4-16B	67-4.6-4-16C	Rayon silk or acetate silk



Caution:

1. This product is for single use only, discard after single use.
2. Do not use it if the package be affected with damp, or the product went moldy.
3. Do not use after expiry date. If over the expiry date, please discard it.
4. Single use only. Latex free.

Storage:

Store in a cool, dry and dark place. Avoid direct sunlight.

Silk Tape



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



EC REP Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com