


# Adhesive Plaster




## **Instructions For Use:**


- Take out the adhesive plaster from the package.
- Peel off the back release paper, exposing the adhesive surface.
- The adhesive surface paste on the wound.
- Used for all kinds of acute and chronic small wound.

## **Features:**


- Material: flexible fabric with absorbent pad, cotton fabric with absorbent pad, non-woven with absorbent pad, PU/PE/PVC with absorbent pad.
- Size:  $\Phi$ 22mm,  $\Phi$ 25mm, 40×10mm, 38x19mm, 55x19mm, 60x19mm, 65×19mm, 70x19mm, 72x19mm, 76x19mm, 76x25mm.
- Non-stick pad for quicker healing, non-adherent to wounds.
- Good breathable and anti-bacterium.




Size	Ref	Material
Φ 22mm	82-F-1-16A	Flexible fabric
Φ 25mm	82-F-2-16A	Flexible fabric
40mm×10mm	82-F-3-16A	Flexible fabric
38mm×19mm	82-F-4-16A	Flexible fabric
55mm×19mm	82-F-5-16A	Flexible fabric
60mm×19mm	82-F-6-16A	Flexible fabric
65mm×19mm	82-F-7-16A	Flexible fabric
70mm×19mm	82-F-8-16A	Flexible fabric
72mm×19mm	82-F-9-16A	Flexible fabric
76mm×19mm	82-F-10-16A	Flexible fabric
76mm×25mm	82-F-11-16A	Flexible fabric



Size	Ref	Material
Φ 22mm	82-C-1-16A	Cotton fabric
Φ 25mm	82-C-2-16A	Cotton fabric
40mm×10mm	82-C-3-16A	Cotton fabric
38mm×19mm	82-C-4-16A	Cotton fabric
55mm×19mm	82-C-5-16A	Cotton fabric
60mm×19mm	82-C-6-16A	Cotton fabric
65mm×19mm	82-C-7-16A	Cotton fabric
70mm×19mm	82-C-8-16A	Cotton fabric
72mm×19mm	82-C-9-16A	Cotton fabric
76mm×19mm	82-C-10-16A	Cotton fabric
76mm×25mm	82-C-11-16A	Cotton fabric



Size	Ref	Material
Φ 22mm	82-N-1-16A	Non-woven
Φ 25mm	82-N-2-16A	Non-woven
40mm×10mm	82-N-3-16A	Non-woven
38mm×19mm	82-N-4-16A	Non-woven
55mmx19mm	82-N-5-16A	Non-woven
60mm×19mm	82-N-6-16A	Non-woven
65mm×19mm	82-N-7-16A	Non-woven
70mm×19mm	82-N-8-16A	Non-woven
72mm×19mm	82-N-9-16A	Non-woven
76mm×19mm	82-N-10-16A	Non-woven
76mmx25mm	82-N-11-16A	Non-woven



Size	Ref	Material
Φ 22mm	82-P-1-16A	PU/PE/PVC
Φ 25mm	82-P-2-16A	PU/PE/PVC
40mm×10mm	82-P-3-16A	PU/PE/PVC
38mm×19mm	82-P-4-16A	PU/PE/PVC
55mmx19mm	82-P-5-16A	PU/PE/PVC
60mm×19mm	82-P-6-16A	PU/PE/PVC
65mm×19mm	82-P-7-16A	PU/PE/PVC
70mm×19mm	82-P-8-16A	PU/PE/PVC
72mm×19mm	82-P-9-16A	PU/PE/PVC
76mm×19mm	82-P-10-16A	PU/PE/PVC
76mmx25mm	82-P-11-16A	PU/PE/PVC

**Caution:**

1. This product is sterilized for single use only, discard after single use.

Please do not repeat sterilization, do not use it once again.

2. Should be used the product immediately after opening the individual package, please destroy when it was used.

3. Do not use it if the inner package has been previously opened or damaged, or package be affected with damp, or the product went moldy.

4. Do not use after expiry date. If over the expiry date, please discard it.

5. Sterile by ETO. Single use only. Latex free.

**Storage:**

Store in a cool, dry and dark place. Avoid direct sunlight.

# Adhesive Plaster



**Secured Medical Direction UK Co., Ltd.**

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



**EC REP** Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

[www.smd-medical.com](http://www.smd-medical.com)